

# PORTABLE OUTLET



*Plug in power anywhere*

## Simultaneously power your CPAP, Phone & Tablet

### MODEL PO-155W

- 155 Watt Battery
- AC Output: 110v AC 60hz/100w
- Two USB Ports: 5 Volt/2.1 amps
- LED Power Indication Meter
- Meets TSA and FAA standards
- Over 500 charge cycles
- Recharges in less than 4 hours
- Size: 5" wide x 7" long x 1.5" tall
- Weight: 2 lbs.
- 1 Year Warranty



***One battery powers EVERY brand and model of CPAP***

### **No need to stock 12 Volt and 24 Volt Models AGAIN!**

One model, without adapters, powers BMC, DeVilbiss, Fisher & Paykel, Resmed, Respiroics, Transcend & more!

### **The longest lasting, most versatile CPAP Battery**

#### Resmed AirMini AutoSet

14 hours set at CPAP 10cmH20 • 10 hours set at CPAP 15cmH20

#### Resmed Airsense 10 AutoSet

8.5 hours CPAP 10cm • 7.5 hours set at CPAP 15cm

#### Philips Respiroics Dreamstation GO AutoCPAP

8 hours set at CPAP 15cmH20

#### Philips Respiroics Dreamstation AutoCPAP

10 hours set at CPAP 15cmH20 (Humidifier Set at 3)

#### Philips Respiroics System One AutoCPAP

11.5 hours set at CPAP 10cmH20 (Humidifier Set at 3)

#### Somnetics/Transcend Auto Mini

10.5 hours set at CPAP 15cmH20

All units tested on manometer. Results may vary depending on heated humidification, and C-Flex or EPR which drops CPAP pressure on exhalation.



Power the AirMini & 2 USB compatible devices



**1595 Peachtree Pkwy, Suite 204-135 • Cumming, GA 30041**  
**www.PortableOutlet.com • 470-297-5778**